

Date 4/6/2010 11:02:00 AM
Description ECG Serial to USB Adapters
Manager ADMIN

Hot Tip Text**PC Reboots when streaming ECG via USB**

The Universal ECG can connect to a PC through a standard DB9 Serial Port (power boost cable recommended) or a USB Port. The Universal ECG ships standard with a Serial-to-USB adapter.

Issue:

The ATEN Serial-to-USB adapter (pn: 1030-3152) causes the PC to reboot in the middle of the test on 64-bit operating systems.

Solution:

Use the FTDI Serial-to-USB adapter (pn: 1030-3155). Contact QRS Technical Support to purchase.

QRS Technical Support
support@qrsdiagnostic
800-465-8408, ext 945