

Hot Tips Report

ID	0021	Description	TN 04.05 - Restoring an IDMS DB on MSDE 2000
Date	4/26/2005 10:38:01 AM		
Manager	ADMIN	Status	Active

Text **Restoring a Small Office Ed. Database**

The information in this Tech Note applies to:

- Office Medic IDMS Small Office Edition 4.01 and higher.

This technical note is only to be used by the Database Administrator (DBA).

For the purposes of this document the term Small Office Edition refers to Office Medic™ (IDMS) Small Office Edition.

Summary

This technical note provides the steps for reattaching the Office Medic database if MSDE 2000 has been un-installed from the server.

Instructions for Restoring Database

A software utility is supplied to reattach the database in this situation. Perform the following steps:

1. Re-install MSDE 2000 onto the server using the Small Office Edition CD-ROM or download files from the QRS website.
 2. Once MSDE has been re-installed and the server has been re-booted, browse the server and find the MSDEEnterpriseManager.exe file.
 3. Double click on MSDEEnterpriseManager.exe and follow the prompts for reattaching the database.
 4. A "MSDE Database Successfully Attached" message is displayed when the database has been successfully reattached.
-